

**AMR-TV**

introductory slide deck

as of March 2025

# Outline

Installation	3
Launching the application	4
Input data overview	6
Uploading KPC3 ecl719 data	7
• Free-zoom view overview	10
Uploading multi-MLST KPC3 data	12
• Zoomed-out view overview	14
Generating config files	15

# Installation

## Requirements:

- Basic familiarity with a command line interface
- Basic familiarity with Github
- Docker

## Other than that, two simple steps:

1) Clone the repo and adaptagrams submodule:

```
$ git clone git@github.com:ivansg44/AMR-TV.git amr-tv --recurse-submodules
```

2) Build the docker image:

```
$ docker-compose build (may take a long time)
```

# Launching the application

Pretty simple:

```
$ docker-compose up
```

AMR-TV will be available at <http://0.0.0.0:8050/>




(see next slide for what launch page should like)

Upload data


Create config file

Free-zoom

Zoomed out



# Input data overview




Example data files are available in the `kpc3_data/` folder of the AMR-TV repo

- These files describe KPC bacteria sampled across British Columbia, Canada


Inputted together:


- `kpc3_data.csv`, `kpc3_config.json`
  - Large dataset of KPC bacteria across multiple MLST
- `kpc3_{MLST}_data.csv`, `kpc3_{MLST}_config.json`, `kpc3_{MLST}_snp_matrix.csv`
  - Several smaller datasets for KPC bacteria with a specific MLST


Uploading KPC3 ecl719 data











primary y-axis  
attribute:

Submitter

secondary y-axis  
attributes:

PrimaryID

Left\_flanks;Right\_flanks

Relaxes

Incompatibility\_Groups

Surname

Environmental

Clinical

Click to filter out  
nodes by shape

potential hgt

Click to filter out  
specific link types

potential clonal spread

Slider to filter potential  
clonal spread links by  
custom cgSNP values

Weight=3 6

Predicted\_mobility

conjugative

Click to filter out  
nodes by color


Note:


node shapes == environmental vs  
clinical samples

node colors == conjugative vs  
mobilizable vs non-mobilizable  
plasmids

In this dataset, we only have  
conjugative plasmids

Uploading multi-MLST KPC3 data





# Generating config files

Standardized JSON vocabulary

You can use the “create config file” modal to help generate a config file

- Requires an existing tabular dataset to generate the correct vocabulary

Multiple “Help” buttons to guide you through the process

### Create config file

kpc3\_data.csv

Comma

Help

**Sample ID field** is the field from your data that contains your sample identifiers. Hovering over a node will display the Sample ID value corresponding to that node. If you upload a matrix, the row and column headers must correspond to this Sample ID field.

Sample ID field:

SampleID

Help

Date field:

Date input format:

Generate config file